Accessories for NIRS™ DA1650











Accessories for NIRS™ DA1650

Large range of accessories available for NIRSTM DA1650. Easy to use and enabling accurate and fast analyses of both solid and liquid samples. Ideal to help you manage your daily analysis from raw material intake to final product control.

	Accessory		Order no.
O	Large cup for pellets, whole grains or other inhomogeneous samples	Single cup	60056582
		Set of 5 cups	60061154
	Small cup for milled or homogeneous samples. Use with holder for small sample cup	Single cup	60053059
		Set of 10 cups	60053921
	Holder for small sample cup		60083686
9	Slurry cup for liquids and pastes. Use with a gold reflector	Single cup	60053173
		Set of 5 cup	60071725
	6 mm gold reflector for liquid oil analyses. Use with slurry cup		60060546
	Stainless steel lids for covering cups	Large / extended cup	10014368
		Small cup	10014367
	Large cup extension for forage or other bulky samples		60063299
00	Adaptor for petri dishes	89 mm petri dishes	60047398
		58 mm petri dishes	60047410